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***Measurement and instrumentation***

**Ref :  
293 057**

**Français – p 1**

**English – p 17**

**Version : 9001**

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**NOTES**


- ① On – off switch
- ② Power-on indicator lamp
- ③ Voltage selector switch:  $3 \times 6 \text{ V} - 3 \times 9 \text{ V}$  RMS between the live and neutral
- ④ Frequency selector switch: 0.1 Hz - 1 Hz - 50 Hz -150 Hz and “0.1 Hz” with manual displacement
- ⑤ Output terminals

## 1 Measuring the phase to ground and phase to phase voltages

In a three-phase installation there are three lives and one neutral wire.

The voltage between the phase and the neutral is called phase to ground voltage; its effective value is noted as:  $V$ .

The voltage between the two phases is called the line voltage; its effective value is noted as:  $U$ .

## 1.1 Objective

To check the relation between phase to ground voltage  $V$  and line voltage  $U$ :

$$U = V \times \sqrt{3} .$$

## 1.2 Equipment

- An extra-low voltage three-phase generator;
- Three digital voltmeters
- Six connecting wires.

## 1.3 Measurements

- Measure the phase to ground and line voltages; complete the table of measurements.

Example of results:

Rated capacity	Phase to ground voltages			Line voltages		
	$V_{1N}$	$V_{2N}$	$V_{3N}$	$U_{12}$	$U_{23}$	$U_{31}$
$3 \times 6 \text{ V}$						
$3 \times 9 \text{ V}$	$V_{1N} = 8.9 \text{ V}$	$V_{2N} = 8.9 \text{ V}$	$V_{3N} = 8.9 \text{ V}$	$U_{12} = 15.5 \text{ V}$	$U_{23} = 15.3 \text{ V}$	$U_{31} = 15.4 \text{ V}$

- Compare for each of the two rated capacities the values of the phase to ground voltages  $V$  then the line voltages  $U$ .
- Check the relation  $U = V \times \sqrt{3}$

## 2 Viewing the voltages using ESAO4 and Générís

### 2.1 Objective


- To view in three-phase sinusoidal mode: The voltages between the phase and the neutral  
The voltages between the phases.
- To verify that a three-phase voltage system is a system with three sinusoidal voltages having the same amplitude, same frequency and offset one from the other by a third of a period
- To verify that the line voltages are leading by  $\frac{\pi}{6}$  from the phase to ground voltages
- To construct Fresnel's diagram for the phase to ground and line voltages (see 4).

## 2.2 Equipment

- An extra-low voltage three-phase generator (9V - 1 Hz)
- A computer with Générís installed
- ESAO4 console
- Three voltmeter adaptors
- Six connecting wires.

## Phase to ground voltages


## 2.3 Set-up

## 2.4 Experiment

- Perform the assembly.
- Start the computer; click on the Générís icon (the indicator lamp on the console comes on).
- Select the Test and calibrate icon
  - Select: Channel: **Voltmeter /**  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
  - Select: Channel: **Voltmeter //**  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
  - Select: Channel: **Voltmeter ///**  
Rated capacity:  25 V  
 Instantaneous

- Adjust the zero
- Click on OK
- Select the Parameters icon
  - Validate: Acquisition mode: f(t) ▼  
 Acquisition time: 2 s ▼  
 Number of points: 501 ▼  
 Synchronisation
  - Select: Channels on the ordinate: Voltmeter/  
 Voltmeter//  
 Voltmeter///
  - Click on Synchrono...
  - Select: Channel: Voltmeter/ ▼  
 Level: 0  
 Synchronisation  Increasing
  - Click on OK
  - Click on OK
- Select the Launch icon
  - Select: Magnitude: v1  
 Names of the curves to be acquired: v1(t)  
 Magnitude: v2  
 Names of the curves to be acquired: v2(t)  
 Magnitude: v3  
 Names of the curves to be acquired: v3(t)
  - Click on Acquisition

Carry out the same search for:

$$v_2(t) = a \cdot \sin(\omega t + b)$$


and

$$v_3(t) = a \cdot \sin(\omega t + b)$$

We obtain for example:

Parameters a  12.8  
 w  6.29  
 b  2.02  
 Standard deviation S =  
 201mV

Parameters a  12.9  
 w  6.27  
 b  4.14  
 Standard deviation S =  
 209mV

Check that  $v_1(t) + v_2(t) + v_3(t) = 0$  ( $\vec{V}_1 + \vec{V}_2 + \vec{V}_3 = \vec{0}$ ).

- Select Processing then Formula
  - Select: Curve to be processed: v1(t) ▼
 

Ordinate:	Name: v
	S.I. Units: V
Formula	v = : v1+v2+v3
Resulting curve:	v(t)
  - Click on Plot then Save and then Close
- Select Processing then Regression
  - Select: Curve to be processed: v(t) ▼
 


Processing interval	Start: 0
	End: 2
Equation	V(t)
  - Click on Plot

We read for example:

$$23.9E-3 t + 182E-3$$


Correlation coefficient:

0.2%

## Line voltages

### 2.6 Set-up


### 2.7 Experiment

- Perform the assembly.
- Select the Test and calibrate icon
  - Select: Channel: Voltmeter /  
 Rated capacity: Ⓞ 25 V  
 Ⓞ Instantaneous

- Adjust the zero
- Click on OK
- Select: Channel: Voltmeter//  
Rated capacity:  25 V  
 Instantaneous
- Adjust the zero
- Click on OK
- Select: Channel: Voltmeter///  
Rated capacity:  25 V  
 Instantaneous
- Adjust the zero
- Click on OK
- Select the Parameters icon
  - Validate: Acquisition mode, f(t) ▼  
Acquisition time: 2 s ▼  
Number of points: 501 ▼  
 Synchronisation
  - Select: Channels on the ordinate: Voltmeter/  
Voltmeter//  
Voltmeter///
  - Click on Synchro...
  - Select: Channel: Voltmeter/ ▼  
Level: 0  
Synchronisation   
 Increasing
  - Click on OK
  - Click on OK
- Select the Launch icon
  - Select: Magnitude: u12  
Names of the curves to be acquired: u12(t)  
Magnitude: u23  
Names of the curves to be acquired: u23(t)  
Magnitude: u31  
Names of the curves to be acquired: u31(t)

- Click on Acquisition

## 2.8 Processing the curves


- Click on the Representation of curves icon
  - Select: Scales:  Automatic.
  - Click on the Colours & Styles icon to modify the colour and the lines of the curves.

Represent  $u_{12}(\pi)$  ,  $u_{23}(\pi)$  and  $u_{31}(\pi)$ .

- Select Processing then Mapping
  - Select: Curve to be mapped: u12(t) ▼
 

Abscissa	S.I. Units:	pi
	Expression:	2*t
Ordinate:	Name:	u12
	S.I. Units:	V
	Expression	u12
Mapped curve:		u12(phi)
  - Click on Plot then Save and then Close

Do the same with the other curves  $u_{23}(t)$  and  $u_{31}(t)$ .


- Click on the Pointer icon
  - Use the cross hair to measure:
    - The period of each of the line voltages
    - The offset from one line voltage to the other
    - The maximum value of each of the line voltages.

Use Zoom for a more accurate reading.

## Phase to ground and line voltages

### 2.9 Set-up


### 2.10 Experiment

- Perform the assembly.
- Start the computer; click on the Générés icon (the indicator lamp on the console comes on).
- Select the Test and calibrate icon
  - Select: Channel: Voltmeter /  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
  - Select: Channel: Voltmeter//  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
  - Select: Channel: Voltmeter///  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
- Select the Parameters icon
  - Validate      Acquisition mode: f(t) ▼  
Acquisition time: 2 ms ▼  
Number of points: 501 ▼  
 Synchronisation
  - Select: Channels on the ordinate: Voltmeter/  
Voltmeter//  
Voltmeter///


- Click on Synchro...
- Select: Channel: **Voltmeter/** ▼  
Level: 0  
Synchronisation      Increasing
- Click on OK
- Click on OK
- Select the Launch icon
  - Select: Magnitude: v1  
Names of the curves to be acquired: v1(t)  
Magnitude: v2  
Names of the curves to be acquired: v2(t)  
Magnitude: u12  
Names of the curves to be acquired: u12(t)
  - Click on Acquisition

## 2.11 Processing the curves

- Click on the Representation of curves icon
  - Select: Scales:  Automatic.
- Click on the Colours & Styles icon to modify the colour and the lines of the curves.
- Click on the Pointer icon
  - Use the cross hair to measure:
    - The period for each of the phase to ground and line voltages
    - The phase shift in line voltages  $u_{12}$  in relation to phase to ground voltage  $v_1$
    - The maximum value of each of the voltages.

Use Zoom for a more accurate reading.

- Select Processing then Formula
  - Select: Curve to be processed: v1(t) ▼  
 Ordinate: Name: U12  
 S.I. Units: V  
 Formula U = : v1-v2  
 Resulting curve: U12(t)

- Click on Plot then on Save

Compare  $u_{12}(t)$  and  $U_{12}(t)$ .

### 3 Viewing the voltages on an oscilloscope

#### 3.1 Objective (see 2)


#### 3.2 Equipment

- An extra-low voltage three-phase generator
- A dual-trace oscilloscope
- Four connecting wires.

#### 3.3 Experiment


Work with a 150 Hz frequency to obtain a stable signal.

- ❖ View  $v_1$  and  $v_2$

- Press the INVERT Y2 button to view  $-v_2$
- Press the ADD button to view  $v_1 - v_2$
- Plot the oscillogram

- ❖ View  $u_{12}$

Compare the two signals  $v_1 - v_2$  and  $u_{12}$ .

To view the three voltages, use the two-way switch for phases 2 and 3.

## 4 Fresnel representation of voltages

We associate to every sinusoidal signal  $u(t)$  a Fresnel vector  $\vec{U}$  that rotates at an angular velocity  $\omega$ .


Represent the Fresnel vectors  $\vec{V}_1$ ,  $\vec{V}_2$  and  $\vec{V}_3$  associated to the three phase to ground voltages  $v_{1N}$ ,  $v_{2N}$  and  $v_{3N}$ .

Represent the Fresnel vector  $\vec{U}_{12}$  associated to the line voltage  $u_{12}$  knowing that:  $u_{12} = v_{1N} + v_{N2}$

That is  $\vec{U}_{12} = \vec{V}_1 - \vec{V}_2$

Graphically check the relation between the phase to ground and line voltage.

Such a system is called balanced.

## 5 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries